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PRODUCENT:

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MANUFACTURER:

Oxylne Sp. z o.o. Poland, 95-200 Pabianice, ul. Pilsudskiego 23.

Manufactured according to PPE Regulation (EU) 2016/425, EN 149:2001+A1:2009.

QUALITY IN ACCORDANCE WITH STANDARDS AND REGULATIONS:

• **EN14683:2019 + AC:2019 • Regulation (EU) 2017/745** – CLASS I TYPE IIR MEDICAL DEVICE. Before reading the following information, please check that the product is a medical device. Information is provided on the packaging. **For medical devices**, this instruction applies to models XR250 FFP2 NR, XR350 FFP3 NR, X310 FFP3 NR, X210 FFP2 NRD.

• **EN149:2001+A1:2009 • Regulation (UE) 2016/425** – PERSONAL PROTECTIVE EQUIPMENT CATEGORY III. According to EN 149:2001+A1:2009 standard disposable respirators against solid and liquid aerosols are classified according to one of three categories: FFP1, FFP2, and FFP3. Before reading the following information, please check the category to which the respirator belongs - this is indicated on the packaging and on the respirator.

MARKING OF THE HALF MASK: The marking printed on each respirator includes the following: Example marking: **XR* 250** - respirator symbol, **V** - exhalation valve, **FFP1, FFP2, FFP3** - filtration class, **NR** - single use respirator - can be used for a maximum of one work shift, **R** - reusable respirator, **D** - this respirator has passed the optional dolomite dust clogging test for more free and comfortable breathing; **EN 149:2001+A1:2009** - the standard designation for PPE, **EN14683:2019+AC:2019** - standard designation for medical devices, **EXP YYY.YY.DD.MM.XXXXXX** - expiry date and sequence number forms a unique product number.

* **X** - Half Masks Cup Shaped, **XF** - flat half mask, **XR** - folding flat half mask, **XS** - half mask with foam sealing, **XCS** - half mask with foam sealing and active carbon layer, **XC** - half mask with active carbon layer.

MEANING OF PICTOGRAMS: 1. product code 2. See user information 3. Single-use product, do not reuse 4. non-sterile product 5. Maximum relative humidity level at storage 6. storage temperature range 7. expiry date 8. Serial number 9. CE mark in accordance with the Regulation of the European Parliament and of the Council (EU) 2016/425 of March 9, 2016 and Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 and 4-digit number of the notified body performing surveillance according to module C2 10. medical device

PICTOGRAMS: 1. 2. 3. 4. 5. 6. 7. 8. 9. 10.

APPLICATIONS: 1. **MEDICAL DEVICE:** The half-mask is a non-sterile disposable medical device, type II R, complying with EN 14683:2019 with a limited time of use of maximum 8 hours, designed for protection in locations exposed to contact with infectious material. Bacteria filtration efficiency $BFE \geq 98.0\%$. The half-mask is designed to cover mouth, nose and chin. The product provides a barrier restricting transmission of infectious agents from and to human respiratory system. It prevents spreading of air-borne micro-organisms, reducing a risk of infections by healthcare professionals in order to provide sterile treatment conditions and protect patient from viruses/bacteria, e.g. during medical procedures. The half-mask is designed to be used by healthcare professionals and patients in case of exposure to blood and other potentially infectious body fluids, and emissions of viruses and bacteria to the environment, but is not considered a sufficient barrier for user against infection with virus from the environment. The product can be used in medical environment, where there is a risk of spraying of body fluids 2. **PPE:** Filtering half-mask class: • **FFP1** is designed to protect the respiratory system against particles in solid as well as liquid phase - dust, smoke, mists, as long as the concentration of the dispersed phase does not exceed **4xOEL*** (**APF**=4, NPF=4**), **FFP2** is designed to protect the respiratory

system against particles in solid as well as liquid phase - dust, smoke, mists, as long as the concentration of the dispersed phase does not exceed **10xOEL*** (**APF**=10, NPF=12**), • **FFP3** is designed to protect respiratory system against particles in solid and liquid phase - dust, smoke, mists, provided that the concentration of the dispersed phase does not exceed **30xOEL*** (**APF**=30, NPF=50**). A respirator is a complete piece of protective equipment ready for immediate use. It is individual equipment and must not be shared with other persons. Note: * OEL = Occupational Exposure Limit, APF = Assigner Protection Factor, NPF = Nominal Protection Factor, ** depending on local and national regulations.

CONTRAINDICATIONS: Use with caution if allergic to polypropylene. If hypersensitivity develops, discontinue use. Consult a doctor if necessary. The respirator does not provide oxygen. Do not use in oxygen deficient conditions (less than 17%). Do not use in confined spaces, especially unventilated spaces such as sewers, manholes, reservoirs etc. The respirator does not protect against pollutants in the form of vapours and gases of substances harmful to health and dangerous to life. Do not use if the type, concentration and properties of harmful substances are not known. Do not use during firefighting. The respirator will not provide a seal if worn on an unshaven or bearded face. Do not use in explosive atmosphere.

HOW TO USE THE MASK: Before using the respirator, check its shelf life and technical condition: that its components are not damaged. Damaged respirators and those whose expiry date has expired must not be used. Do not fold or bend the respirator.

FITTING INSTRUCTIONS: 1. Disinfect your hands. 2. Before putting on the respirator, pre-shape the nose clip by tightening it, which improves the subsequent fit of the respirator to the face. 3. Position the facepiece so that it covers the mouth and nose. 4. Put on the retaining straps - for respirators series: X, XC, XS, XCS; so that the lower strap goes around the neck below the ear and the upper strap goes around the back of the head above the ear. For half masks fitted with a sliding head strap, adjust the length of the lower and upper head strap to individual needs. • for half masks of the XR series: over the ears. 5. Check that it is fitted correctly. Put your hands together and hold the canopy. Blow out air firmly. In case of leakage, adjust the position of the canopy, nose clip or ear or head harness straps.

Half mask: X, XC, XS, XCS series

Half mask XF series

Half mask XR series

REMOVAL:

1. Go to dirty zone. 2. When using additional personal protective equipment, e.g. overalls, protective goggles, glasses, etc. take them off in accordance with instructions for taking off a particular product 3. Remove outer gloves and discard them to waste bin 4. For X, XC, XS, XCS and XF Series masks, lift the lower head strap up around your neck, then pick up the upper head strap and pull both straps away from your head. For XR series masks, pull the side tapes aside so that they do not touch the skin. Carefully remove the respirator without allowing the outer part of

the respirator to come into contact with the body. 5. Discard half-mask to waste bin. 6. Take off protective gloves and discard them to waste bin. 7. Disinfect hands using proper disinfection agents. ***NOTES:** If breathing becomes difficult, the mask is damaged or a proper fit cannot be maintained, replace the mask immediately. Careful adherence to the above instructions is an important part of safe mask use. If there is a danger of oxygen deficiency or too high a concentration of a hazardous substance, use equipment that isolates the user from the environment and leave the hazardous area. The respirator is designed to be used for one work shift (max. 8 hours) - applies to respirators marked with NR.

TRANSPORT, STORAGE AND MAINTENANCE: The respirator should be transported and stored in its original packaging, in dry areas and away from light and heat sources. As disposable equipment the respirator cannot be reconditioned or repaired. It must not be washed or cleaned. The respirator should be stored in closed, undamaged manufacturer's packaging in rooms with relative humidity below 70% and at temperatures from -30°C to +50°C. If stored properly the shelf life is 60 months. Protect against aggressive chemicals, moisture and dirt. If the mask has not been contaminated, it can be disposed of in accordance with local waste management regulations. If contaminated, it must be disposed of in accordance with applicable laws and regulations.

WARNING:

1. Failure to follow all instructions and limitations on the use of this product could reduce the effectiveness of respirator and result in sickness or death. 2. Before use, make sure that the respirator has no defects: tears, rips, deformation of the valve membrane, visible dirt or any other damage. If found, the unit is not suitable for use and must be disposed of. 3. It is recommended to check that the product is suitable for the intended use. 4. A properly selected respirator is essential to protect your health. Before using this respirator, consult an Industrial Hygienist or Occupational Safety Professional to determine the suitability for your intended use. 5. This product does not supply oxygen. Use only in adequately ventilated areas containing sufficient oxygen to support life. Do not use this respirator when oxygen concentration is less than 17%. 6. Do not use when concentrations of contaminants are immediately dangerous to health or life. 7. Do not use in explosive atmospheres. 8. Leave work area immediately and return to fresh air if (a) breathing becomes difficult, or (b) dizziness or other distress occurs. 9. Facial hairs or beards and certain facial characteristics may reduce the effectiveness of this respirator. 10. Never alter or modify this respirator in any way. 11. A NR marked respirator may be used **(a) as a single use medical device** or **(b)** for a maximum of one shift (8 hours) as PPE. The R-marked respirator may be used for longer than one shift provided that the inner surface of the respirator is disinfected with a 70% ethanol solution (by spraying or wiping) or exposed to a UV lamp only as PPE. No maintenance is required. After use, the mask should be disposed of in accordance with applicable laws and regulations. 12. Keep respirators in the display box away from direct sunlight until use. 13. Respirators XC series, with a carbon layer, reduce odorous odours. They DO NOT PROTECT against gases with concentrations above 1 OEL. In this case, use a gas mask.

CERTIFICATION AND SUPERVISION **:

Certification according to module B: 1. CIOP-BiP, Czerniakowska 16, 00-701 Warszawa, Poland, Notified Body No. 1437. 2. ITT Certex, Górnicka 30/36, 91-765 Łódź, Notified Body No. 2534. 3. APAVE Sudeurope SAS Centre d'Essais et de certification de Fontaine 17, Boulevard Paul Langevin 38000 FONTAINE - France, Notified Body No 0082

Supervision according to module C2: 1. CIOP-BiP, Czerniakowska 16, 00-701 Warszawa, Poland, Notified Body No 1437. 2. ITT Certex, Górnicka 30/36, 91-765 Łódź, Notified Body No 2534.

* **Notes:** Depending on the model, information on certification and supervision is included in the declaration of conformity.

Declaration of Conformity download at: www.oxylne.eu

PRODUCER/DISTRIBUTOR: details on packaging