
(PL) [Karta katalogowa](#)

(EN) [Data sheet](#)

(DE) [Datenblatt](#)

(FR) [Fiche d'information](#)

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Oxyline Sp. z o.o.
95-200 Pabianice, Poland
st. Piłsudskiego 23
tel.: +48 42 215 10 68 e-
mail: oxyline@oxyline.eu


EUROPEAN
PRODUCT

MADE IN
POLAND

OXY[®] XCS 320 V FFP3 NR D

Certyfikat: UE/832/2022/1437, wydanie 1

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
Certificate: UE/832/2022/1437, edition 1

Date of issue: 19 December 2022

The **XCS 320 V FFP3 NR D** filtering half mask is designed to protect the respiratory system against particulate aerosols, water-based aerosols (dust, smoke) and aerosols with liquid phase dispersion (mists) if the concentration of the dispersed phase does not exceed 30 x OEL.

DESCRIPTION:

The filtering half mask **XCS 320 V FFP3 NR D** is composed of the following elements :

- A multi-layered filtration material: polypropylene;
- Inner layer of activated carbon;
- A nose clip to shape the half mask at the nose;
- Exhalation valve;
- Header tapes made of braided rubber thread with plastic fasteners;
- Plastic fastenings of the head bands;
- The internal layer made of polyurethane, to improve the seal and fit to the face.

The half mask is designed in such a way as to enable easy breathing throughout the work shift. The anatomical shape and the nose clip, as well as the internal sealing foam, make the half mask easy to fit to most face shapes, so that the necessary tightness can be ensured.

The half mask **XCS 320 V FFP3 NR D** is designed to protect respiratory system against harmful effect of dust, solid and liquid aerosols when OEL is $\leq 0,05 \text{ mg/m}^3$ and the concentration of dispersed phase of aerosol does not exceed 30 x OEL (Occupational Exposure Limit), 30 x APF (Assigned Protection Factor), 50 x NPF (Nominal Protection Factor). The activated carbon layer reduces the negative effects associated with odours or fume emissions and can be used in environments where there are organic or inorganic vapours and gases with concentrations below 1 x OEL.

EXAMPLES OF APPLICATION

high concentrations of respirable dusts, welding and soldering, waste sorting, protection against dusts containing beryllium, antimony, arsenic, cadmium, cobalt, nickel, radium, strychnine

HOW IT WORKS

The filtering half mask is composed mostly of the face part made of filtering material and accessories such as headbands, or exhalation valve, depending on the model. When air is drawn in, it passes through the filtration material where it is cleansed before being inhaled. Exhaled air passes through filtration material (in the masks without a valve) or through both the exhalation valve and the filtration material (in models with a valve). The cup of the mask should be well adjusted to the user's face.

Respiratory protection

REQUIREMENTS:

OXYLINE half masks comply with the following:

- harmonised European standard EN 149:2001+A1: 2009 "Respiratory protective devices - Filtering half masks to protect against particles. Requirements, testing, marking";
- in conformity with the relevant Union harmonisation legislation: Regulation (EU) 2016/425 of the European Parliament and of the Council of 9 March 2016 on personal protective equipment and repealing Council Directive 89/686/EEC (Text with EEA relevance).

CONTRAINDICATIONS

It does not ensure protection of the respiratory system if there is a lack of oxygen (below 17%). It should not be used in spaces with limited cubic volume, in particular non-ventilated spaces, such as sewers, wells, tanks, etc. The half mask does not provide protection against pollution in the form of gas fumes or mists of substances that are harmful to human health and hazardous to life. Do not use the half mask if the type, characteristics and concentration of the harmful substances are unknown. Do not use the half mask when extinguishing fires. The half mask does not ensure tightness if worn on an unshaven or bearded face.


FUNCTIONAL PARAMETERS OF THE HALF MASK

Class (according to EN 149:2001+A1: 2009)	FFP3
Penetration of filtering material by sodium chloride aerosol or oil mist	≤ 1%
Total leakage	≤ 2%
Initial inhalation resistance at a flow of 95 l/min	≤ 300 Pa
Initial exhalation resistance 160 l/min	≤ 300 Pa
Breathing resistance at the end of clogging test with dolomite dust, at a flow of	95 l/min (inhalation) ≤ 700 Pa
	160 l/min (exhalation) ≤ 300 Pa

USE AND STORAGE

The half masks should be stored at a temperature of -30°C to +50°C and humidity below 70%. Before the half mask is used, its technical condition should be checked, i.e. whether the elements are not damaged. Damaged or expired half masks must not be used. The half mask should not be folded or bent. A respirator marked NR must not be used for more than one 8-hour working shift. In order to ensure the best possible fit on the face, the half masks should be put on and adjusted in the following manner:

1. Disinfect your hands.
2. Before putting on the half mask, form the nose clip by tightening,
3. Place the mask over the face to cover the mouth and the nose;
4. Put the head bands on in such a way as to make the lower band pass around the nape of the neck below the ear, and the upper band pass around the back of the head above the ear; the length of upper and lower band can be adjusted;
5. Further adjust the nose clip to ensure tightness
6. Check that you have the correct mounting. Press your hands and hold the dome of the mask. Exhale energetically; if there is any looseness adjust the position of the dome, the nose clamp or headbands.
7. Shelf life 60 months. The expiry date on the product.

During the use of the respirator, there is an increase in breathing resistance caused by dust deposition. If the user finds that the resistance has increased significantly, the respirator must be replaced with a new one. The respirator must not be used for more than one 8-hour working shift.

ADJUSTING THE HEAD BANDS:

To adjust the head bands, put the strap between the plastic hooks as shown in illustration 2. Looping it a number of times between the hooks will shorten the band and tighten the mask when it is put on. In order to extend the length again just take the band off the hooks.


For more precise adjustment of the head band, wind or unwind the band on/off each hook, as shown in illustration 3.

Notified body Central Institute for Labor Protection - National Research Institute, ul. Czerniakowska 16, 00-701 Warsaw, - No. 1437.

Oxyline Sp. z o.o.
95-200 Pabianice, Poland
st. Piłsudskiego 23
tel.: +48 42 215 10 68
e-mail: oxyline@oxyline.eu

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