### Ochrona zgodna z normami

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# OxyChem C210

Protection Coverall, Cat. III, Type 5, 6. Medical product Class I

## Certification

|                                                                               |                                                   |                                                       |                          |                                      |                                            |
|-------------------------------------------------------------------------------|---------------------------------------------------|-------------------------------------------------------|--------------------------|--------------------------------------|--------------------------------------------|
| Type 6 - Protective clothing against liquid chemicals - low-pressure splashes | Type 5 - Protective clothing against particulates | Protective clothing against radioactive contamination | Electrostatic properties | Protection against infectious agents | Protective clothing - General requirements |
|                                                                               |                                                   |                                                       |                          |                                      |                                            |
| EN 13034:2005 +A1:2009                                                        | EN ISO 13982-1:2004 +A1:2010                      | EN 1073-2:2002                                        | EN 1149-5:2018           | EN 14126:2003 +AC:2004               | EN ISO 13688:2013                          |

## Lightweight and tough

The **OxyChem C210** protective suit is designed for protection in areas exposed to infectious materials and provides limited type 5/6 anti-chemical protection. It is made of non-woven fabric laminated with a microporous film (**MPFL**) of 63g/m<sup>2</sup> +/- 2g/m<sup>2</sup>. Outer layer: polyethylene film. Inner layer: polypropylene fibres. Non-woven fabric (**MPFL**) is resistant to permeation infectious agents. **OxyChem C210** also provides protection against radioactive contamination and has anti-electrostatic properties. Its construction is designed to provide the wearer with the highest level of safety, ergonomics and comfort during work.

## Properties

- Spunbond fabric laminated with microporous film (**MPFL**), 63g/m<sup>2</sup> +/- 2g/m<sup>2</sup>
- Three-panel hood
- Two-way zip with fabric patch
- Two-piece crotch cut
- Elastic cuffs, legs, waist and hood
- Elastic thumb loop
- Durable overlock seams on the inside

## Application for jobs

- When working with asbestos
- When painting and spray-painting
- When laying insulation and earthworks
- When working on water supply systems, in sewers
- In the pharmaceutical industry
- When working with resin
- Machine maintenance
- During light spraying work, e.g. in agriculture
- In the automotive industry
- Grinding and polishing
- In the pharmaceutical and medical industries

## Packing

The suit is individually wrapped in a sealed foil package.

- Quantity in a carton - 50 pcs.
- Quantity of overalls on a pallet - 1000 pcs.
- Quantity of cartons on a pallet - 20 pcs.
- Gross weight of cardboard - 12,8 kg
- Dimensions of the cardboard - 48 cm x 28 cm x 54 cm

## Technical parameters

Type of protection obtained according to EN 14126 (Protective clothing - Requirements and test methods for protective clothing against infectious agents) - Performance levels.

| Standard                      | Title                                                                                                                                                                                                                  | Retrieved from type protection |
|-------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------|
| EN 13034:2005+A1:2009         | Protective clothing against liquid chemicals - Performance requirements for chemical protective clothing offering limited protective performance against liquid chemicals (Type 6 and Type PB[6] of clothing)          | Type 6B                        |
| EN ISO 13982-1:2004 + A1:2010 | Protective clothing for use against solid particulates - Part 1: Performance requirements for chemical protective clothing providing protection to the full body against airborne solid particulates (Type 5 clothing) | Type 5B                        |

Performed examinations and tests in accordance with EN 14126 (test for protection against infectious agents):

| Standard test                                                                                                                                            | Result obtained | Classification by result |
|----------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|--------------------------|
| Pressure test, resistance to penetration of blood and body fluids with synthetic blood (EN ISO 16603 + EN 14126 par. 4.1.4.1)                            | 20kPa           | Class 6                  |
| Resistance to pathogen penetration from blood using bacteriophage (Fag Fi-X174) (resistance to viral penetration) (EN ISO 16604 + EN 14126 par. 4.1.4.1) | > 75 min        | Class 6                  |
| Resistance to bacterial penetration (EN ISO 22610 + 14126 par. 4.1.4.2)                                                                                  |                 |                          |
| Resistance to penetration of biologically contaminated aerosols (EN ISO 22611 + EN 14126 par. 4.1.4.3)                                                   | Log10 CFU > 5   | Class 3                  |
| Resistance to penetration of contaminated dust (EN ISO 22612 + 14126 par. 4.1.4.4)                                                                       | Log10 CFU < 1   | Class 3                  |

Performed examinations and tests in accordance with EN 1073-2 (protection against contamination with radioactive particles):

| Standard test                                       | Result obtained        | Classification by result |
|-----------------------------------------------------|------------------------|--------------------------|
| Total internal leakage (EN ISO 13982-2 + EN 1073-2) | TILE < 3%<br>TILA < 2% | Class 2                  |

Parameters relating to the strength of the garment (physical properties)

| Standard test                                                                                                           | Result obtained                       | Classification by result |
|-------------------------------------------------------------------------------------------------------------------------|---------------------------------------|--------------------------|
| Determination of resistance to abrasion of materials (EN 530 met. 2 + EN 14325 par. 4.4.1 + EN 13034 par. 4.1)          | > 1000                                | Class 4                  |
| Determination of resistance to flexural fracture (EN ISO 7854 met. B + EN 14325 par. 4.5 + EN ISO 13982-1 par. 4.1)     | > 100 000                             | Class 4                  |
| Determination of trapezoidal tearing strength (EN ISO 9073-4 + EN 14325 par. 4.7 + EN 13034 par. 4.1)                   | Longitudinal > 60<br>Transverse > 40  | Class 4/3                |
| Determination of resistance to push-out (EN ISO 13938-1 + EC 1 + EN 14126 par. 4.1.2)                                   | > 320                                 | Class 4                  |
| Tensile strength and elongation. Strip method (EN ISO 13934-1 + EN 14325 par. 4.9 + 13034 par. 4.1)                     | Longitudinal > 100<br>Transverse > 60 | Class 3/2                |
| Puncture resistance of materials (EN 863 + 1073-2 par. 4.2)                                                             | > 50                                  | Class 3                  |
| Electrostatic properties test Surface resistivity (EN 1149-1 + EN 1149-5 par. 4.2.1)                                    | < 2,5 x 10 <sup>9</sup> Ω             |                          |
| Seam strength (max. breaking load) using the grab method (EN 13935-2 + EN 14325 para. 5.5 + EN ISO 13982-1 para. 4.2.2) | > 100                                 | Class 3                  |
| Flammability tests (ignition and flammability) (13274-4 met 3 + EN 1073-2 par. 4.2)                                     |                                       | Test passed              |

Examinations and tests according to EN 13034 (protection against liquid chemicals)

| Standard test                                                                                               | Result obtained                                                                             |
|-------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
| Test of resistance against chemicals (EN ISO EN ISO 6530 + EN 14325 par. 4.12 and 4.13 + EN 13034 par. 4.1) | H <sub>2</sub> SO <sub>4</sub> (30%) < 1<br>NaOH (10%) < 1<br>o-Ksylan < 1<br>N-butanol < 1 |
| Resistance to penetration by sprayed liquid (Spray test) (EN ISO 17491-4 + EN 14605 par. 4.3.4)             | Test passed                                                                                 |

Additional research:

| Standard test                                                                                                                                                                 | Result obtained                             |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|
| Determination of pH of aqueous extracts (mean PH value of extracts 20 and 30) (EN ISO 3071 + EN ISO 13688 para. 4.2)                                                          | 3,5 > pH > 9,5                              |
| Cytotoxicity tests (Mean % cell survival after incubation) (EN ISO 10993-5)                                                                                                   | 94.47%<br>- is not cytotoxic                |
| Dermatological test (Regulation (EC) No 1223/2009 of the European Parliament and of the Council of 30 November 2009 on cosmetic products) (irritation index X <sub>mn</sub> ) | Non-irritant product<br>X <sub>mn</sub> = 0 |

Assessment of Declared Product Properties - Protective Suit OxyChem C210/C310 (by J.S. Hamilton Poland Sp. z O.O. - No. 99014/21/CGDA of 22.03.2021)

On the basis of the results of the conducted performance test, the properties of the product KOMBINEZON OXYCHEM C210/C310 declared by the manufacturer were evaluated as follows:

|                                                                    |                         |
|--------------------------------------------------------------------|-------------------------|
| The product allows for easy dressing and removal.                  | 100% positive responses |
| The product is convenient to use.                                  | 96% positive responses  |
| The product does not put undue pressure on the skin (welts).       | 92% positive responses  |
| The product does not cause reddening at the contact with the skin. | 92% positive responses  |
| The product provides good thermal comfort while working.           | 100% positive responses |
| The product is airtight.                                           | 100% positive responses |
| The product adheres well at the cuffs, legs, hood.                 | 96% positive responses  |
| The product allows freedom of movement while working.              | 100% positive responses |

In addition, it was assessed:

|                                                   |                         |
|---------------------------------------------------|-------------------------|
| Proper application of the product.                | 100% positive responses |
| Suitable product texture.                         | 96% positive responses  |
| Overall product rating.                           | 84% positive responses  |
| Intention to purchase the product.                | 88% positive responses  |
| Willingness to recommend the product to a friend. | 88% positive responses  |

|                      |                                                                                                                                                                                                                 |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| EN ISO 15223-1:2016  | Medical devices - Symbols to be used on the labelling of medical devices, in conjunction with Medical devices - Symbols to be used for their labelling, marking and information - Part 1: General requirements. |
| EN 1041:2008+A1:2013 | Information supplied by the manufacturer of medical devices.                                                                                                                                                    |
| EN ISO 14971:2019    | Medical devices. Application of risk management to medical devices.                                                                                                                                             |

The EU Declaration of Conformity is available at: <https://www.oxyline.eu/deklaracje-kombinezony.html>

The OxyChem C210 one-piece protective suit is a Class I, non-sterile medical device which complies with the requirements of Regulation 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC (Text with EEA relevance.) and complies with the Medical Devices Act of 20.05.2010. It is also a Category III PPE in accordance with Regulation EU 2016/425 of the European Parliament and of the Council of 9 March 2016 on personal protective equipment and repealing Council Directive 89/686/EEC (Text with EEA relevance). Product designed and placed on the market in accordance with Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 and Regulation EU 2016/425 of the European Parliament and of the Council of 9 March 2016.

| Sizing | S         | M         | L         | XL        | XXL       | XXXL      |
|--------|-----------|-----------|-----------|-----------|-----------|-----------|
| Height | 164 - 170 | 170 - 176 | 176 - 182 | 182 - 188 | 189 - 194 | 194 - 200 |
| Chest  | 84 - 92   | 92 - 100  | 100 - 108 | 108 - 117 | 116 - 124 | 124 - 132 |

Exposure to certain chemicals or to high concentrations may require higher protective properties through the suit's material or construction properties. Such areas may be protected by suit types 1,2,3 or 4.

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