



- Light stretcher (5.5kg), ideal for evacuation in confined spaces.
- Once the victim is installed, the stretcher becomes rigid.
- Can be used horizontally and vertically.
- 2 carrying handles on each longitudinal side.
- Set of straps for a quick and secure attachment of the patient.
- Set of straps for lifting.
- **Maximum load:** 150 Kg.
- Provided with a compact and easy-to-carry transport bag.
- **Total weight with bag and straps:** 8 kg.

Material of the components:

Stretcher: Polypropylène.

Straps: Polyester and Polyamide.

Handles: Polyester.

Size:		Weight:
2.5 m x 0.9 m		5.5kg

KRATOS SAFETY

689, chemin du Buclay
38540 Heyrieux – FRANCE

Tel : +33 (0)4 72 48 78 27

Fax : +33 (0)4 72 48 58 32

www.kratossafety.com / info@kratossafety.com

HS CODE

3901 90 90

CARTON SIZE / PIECES

38 x 38 x 95 cm / 1

Satra Technology Centre, N°0321
Wyndham Way, Telford Way, Kettering
Northamptonshire, NN16 8SD, United Kingdom

ARTICLE 11 GB07/73551

SGS United Kingdom Ltd., N°0120
Unit 202B, Worle Parkway,
Weston-super-Mare, BS22 6WA, United Kingdom